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STD Processing: BP 20-450Hz | Env 6Hz - Solar A.

3.1. Flujo de Operaciones (Pipeline)

1. Filtering (High Pass) - Limpieza Offset

- **Type:** IIR High Pass
- **Frequency:** 20 Hz
- **Order:** 4
- **Phase Lag:** Zero Phase (Bidirectional) ← *Crucial para coincidir con filtfilt.*

2. Filtering (Low Pass) - Limpieza Ruido

- **Type:** IIR Low Pass
- **Frequency:** 450 Hz
- **Order:** 4
- **Phase Lag:** Zero Phase (Bidirectional)

3. Rectification - Valor Absoluto

- **Type:** Full Wave

4. Filtering (Low Pass) - Envolvente Lineal

- **Type:** IIR Low Pass
- **Frequency:** 6 Hz
- **Order:** 4
- **Phase Lag:** Zero Phase (Bidirectional)

EMG Signal Processing

Filtering Protocols, Normalization, and Feature
Extraction
(Python vs Noraxon Ultium)

Author: Alejandro Solar Iglesias

Date: January 2026

Version: 1.0 English

4. Introduction

This document details the digital signal processing (DSP) methodology applied to electromyographic recordings. The objective is to ensure mathematical consistency between automated analysis via Python scripts and visual processing in the proprietary Noraxon myoRESEARCH (MR3) software.

The algorithms described here follow international ISEK standards and SENIAM recommendations.

5. Signal Processing in Python

The processing pipeline has been implemented in the `noraxon_analytics.py` library. Critical stages are described below.

5.1. Spectral Cleaning and Offset Removal

The raw signal contains low-frequency noise (cable movement) and high-frequency noise. A zero-phase **4th-order Band-Pass Butterworth** filter is applied (`filtfilt`).

- **Sampling Frequency (Fs):** 2000 Hz (Ultium Hardware).
- **Low Cutoff:** 20 Hz (DC Offset removal and motion artifacts).
- **High Cutoff:** 450 Hz (Anti-aliasing and electrical noise).

5.2. The Linear Envelope

To analyze "movement intention" and muscle synergy, the oscillatory signal is transformed into a smooth amplitude curve.

1. **Rectification:** Absolute value of the signal ($|x|$).
2. **Smoothing:** Low-Pass Filter at **6 Hz**.

5.3. Robust MVC Calculation and Normalization

To prevent spurious noise spikes ($\leq 10\text{ms}$) from distorting normalization, the absolute maximum is not used. A **500ms Moving Average** algorithm is implemented.

5.4. Feature Extraction

Based on the work of *Oleinikov et al. (2018)*, features are extracted in time and frequency domains to feed classification models.

- **Time Domain:** MAV, RMS, Zero Crossing (ZC), Slope Sign Change (SSC).
- **Frequency Domain:** Mean Frequency (MNF) for fatigue, Wavelet Energy (db4, level 4) for MUAP detection.

6. Signal Processing with Noraxon Software (MR3)

To ensure the software visualization matches the Python code results mathematically, the ”**Signal Processing**” module must be configured following this exact order of operations.

This configuration has been saved in the system as:

STD Processing: BP 20-450Hz | Env 6Hz - Solar A.

6.1. Operations Flow (Pipeline)

1. Filtering (High Pass) - Offset Cleaning

- **Type:** IIR High Pass
- **Frequency:** 20 Hz
- **Order:** 4
- **Phase Lag:** Zero Phase (Bidirectional) $\leftarrow$ *Crucial to match filtfilt.*

2. Filtering (Low Pass) - Noise Cleaning

- **Type:** IIR Low Pass
- **Frequency:** 450 Hz
- **Order:** 4
- **Phase Lag:** Zero Phase (Bidirectional)

3. Rectification - Absolute Value

- **Type:** Full Wave

4. Filtering (Low Pass) - Linear Envelope

- **Type:** IIR Low Pass
- **Frequency:** 6 Hz
- **Order:** 4
- **Phase Lag:** Zero Phase (Bidirectional)